

Product datasheet for **MK300363**

Mir496 Mouse qPCR Template Standard (MI0004589)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0003738
Synonyms:	MI0004589; MIMAT0003738; Mir496 mmu-mir-496; mmu-miR-496
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)