

Product datasheet for **MK300383**

Mir568 Mouse qPCR Template Standard (MI0005517)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0004892

Synonyms: MI0005517; MIMAT0004892; Mir568 mmu-mir-568; mmu-miR-568

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)