

Product datasheet for **MK300389**

Mir653 Mouse qPCR Template Standard (MI0005557)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0004943

Synonyms: MI0005557; MIMAT0004943; Mir653 mmu-mir-653; mmu-miR-653

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)