

Product datasheet for **MK300395**

Mir669a-1 Mouse qPCR Template Standard (MI0004523)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0003477
Synonyms:	MI0004523; MIMAT0003477; Mir669a-1 mmu-mir-669a-1; mmu-miR-669a
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)