

## Product datasheet for **MK300433**

### Mir694 Mouse qPCR Template Standard (MI0004664)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** MIMAT0003474

**Synonyms:** MI0004664; MIMAT0003474; Mir694 mmu-mir-694; mmu-miR-694

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)