

Product datasheet for **MK300444**

Mir707 Mouse qPCR Template Standard (MI0004691)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0003497

Synonyms: MI0004691; MIMAT0003497; Mir707 mmu-mir-707; mmu-miR-707

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)