

## Product datasheet for **MK300501**

### Mir125b-2 Mouse qPCR Template Standard (MI0000152)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Template Standards   |
| Forward Sequence: | TCCCTGAGACCCTAACTTGT  |
| Reverse Sequence: | GAACATGTCTGCGTATCTC   |
| ACCN:             | MIMAT0000136  |
| Synonyms:         | MI0000152; MIMAT0000136; Mir125b-2 mmu-mir-125b-2; mmu-miR-125b-5p  |
| Components:       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control:  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment   |
| Storage:          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)