

Product datasheet for **MK300504**

Mir135a-2 Mouse qPCR Template Standard (MI0000715)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TATGGCTTTTATTCTATGTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000147
Synonyms:	MI0000715; MIMAT0000147; Mir135a-2 mmu-mir-135a-2; mmu-miR-135a
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)