

Product datasheet for **MK300505**

Mir138-1 Mouse qPCR Template Standard (MI0000722)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGCTGGTGTGGAATCAGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000150

Synonyms: MI0000722; MIMAT0000150; Mir138-1 mmu-mir-138-1; mmu-miR-138

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)