

## Product datasheet for **MK300509**

### Mir181b-2 Mouse qPCR Template Standard (MI0000823)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** AACATTCATTGCTGTCGGTGG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0000673

**Synonyms:** MI0000823; MIMAT0000673; Mir181b-2 mmu-mir-181b-2; mmu-miR-181b

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)