

Product datasheet for **MK300511**

Mir196a-2 Mouse qPCR Template Standard (MI0000553)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TAGGTAGTTTCATGTTGTTGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000518
Synonyms:	MI0000553; MIMAT0000518; Mir196a-2 mmu-mir-196a-2; mmu-miR-196a
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)