

Product datasheet for **MK300514**

Mir19b-2 Mouse qPCR Template Standard (MI0000546)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | TGTGCAAATCCATGCAAACCTG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0000513 |
| Synonyms: | MI0000546; MIMAT0000513; Mir19b-2 mmu-mir-19b-2; mmu-miR-19b |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)