

Product datasheet for **MK300516**

Mir297a-2 Mouse qPCR Template Standard (MI0000397)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ATGTATGTGTGCATGTGCAT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000375

Synonyms: MI0000397; MIMAT0000375; Mir297a-2 mmu-mir-297a-2; mmu-miR-297a

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)