

Product datasheet for **MK300521**

Mir30c-2 Mouse qPCR Template Standard (MI0000548)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGTAACATCCTACTCTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000514
Synonyms:	MI0000548; MIMAT0000514; Mir30c-2 mmu-mir-30c-2; mmu-miR-30c
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)