

Product datasheet for **MK300524**

Mir466f-2 Mouse qPCR Template Standard (MI0005508)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TACGTGTGTGTGCATGTGCA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0004881
Synonyms:	MI0005508; MIMAT0004881; Mir466f-2 mmu-mir-466f-2; mmu-miR-466f-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)