

Product datasheet for **MK300548**

MIR1251 Mouse qPCR Template Standard (MI0014021)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACTCTAGCTGCCAAAGGCG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0014824

Synonyms: MI0014021; MIMAT0014824; MIR1251 mmu-mir-1251; mmu-miR-1251-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)