

Product datasheet for **MK300596**

MIR3072 Mouse qPCR Template Standard (MI0014035)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGGGACCCCGAGGGAGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0014852

Synonyms: MI0014035; MIMAT0014852; MIR3072 mmu-mir-3072; mmu-miR-3072-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)