

Product datasheet for **MK300623**

MIR3097 Mouse qPCR Template Standard (MI0014090)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CACAGGTGGGAAGTGTGTGT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0014915

Synonyms: MI0014090; MIMAT0014915; MIR3097 mmu-mir-3097; mmu-miR-3097-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)