

Product datasheet for **MK300644**

MIR344d-3 Mouse qPCR Template Standard (MI0004227)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGTCAGGCTAGTGGTTATACT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0014807

Synonyms: MI0004227; MIMAT0014807; MIR344d-3 mmu-mir-344d-3; mmu-miR-344d-3-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)