

Product datasheet for **MK300658**

MIR3473c Mouse qPCR Template Standard (MI0018015)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCTCTCCAGCCCCATAATA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0020614
Synonyms:	MI0018015; MIMAT0020614; MIR3473c mmu-mir-3473c; mmu-miR-3473c
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)