

Product datasheet for **MK300678**

MIR3963 Mouse qPCR Template Standard (MI0016968)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGTATCCCACTTCTGACAC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019341

Synonyms: MI0016968; MIMAT0019341; MIR3963 mmu-mir-3963; mmu-miR-3963

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)