

Product datasheet for **MK300683**

MIR3968 Mouse qPCR Template Standard (MI0016977)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CGAATCCCACTCCAGACAC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019352

Synonyms: MI0016977; MIMAT0019352; MIR3968 mmu-mir-3968; mmu-miR-3968

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)