

Product datasheet for **MK300714**

MIR5100 Mouse qPCR Template Standard (MI0018008)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TCGAATCCCAGCGGTGCC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0020607

Synonyms: MI0018008; MIMAT0020607; MIR5100 mmu-mir-5100; mmu-miR-5100

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)