

Product datasheet for **MK300722**

MIR5112 Mouse qPCR Template Standard (MI0018021)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TAGCTCAGCGGGAGAGC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0020620

Synonyms: MI0018021; MIMAT0020620; MIR5112 mmu-mir-5112; mmu-miR-5112

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)