

Product datasheet for **MK300732**

MIR5124a Mouse qPCR Template Standard (MI0018035)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: GGTCCAGTGACTAAGAGC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0020634

Synonyms: MI0018035; MIMAT0020634; MIR5124a mmu-mir-5124a; mmu-miR-5124a

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)