

Product datasheet for **MK300755**

MIR5623 Mouse qPCR Template Standard (MI0019191)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGGCAGCTGGCAAGTCAGAAT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0022373

Synonyms: MI0019191; MIMAT0022373; MIR5623 mmu-mir-5623; mmu-miR-5623-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)