

Product datasheet for **MK300781**

MIR6346 Mouse qPCR Template Standard (MI0021874)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TAGCAGGTGTCTGCAATGCC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025089

Synonyms: MI0021874; MIMAT0025089; MIR6346 mmu-mir-6346; mmu-miR-6346

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

[View online »](#)