

Product datasheet for **MK300784**

MIR6349 Mouse qPCR Template Standard (MI0021877)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGG AATGTGGGAGGGGCAT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025092

Synonyms: MI0021877; MIMAT0025092; MIR6349 mmu-mir-6349; mmu-miR-6349

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)