

Product datasheet for **MK300786**

MIR6351 Mouse qPCR Template Standard (MI0021879)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TAGCGTTGCTCAGGGAGCA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025094

Synonyms: MI0021879; MIMAT0025094; MIR6351 mmu-mir-6351; mmu-miR-6351

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)