

## Product datasheet for **MK300788**

### **MIR6353 Mouse qPCR Template Standard (MI0021881)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** TAGCAGCACGTATTTATCCT

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0025096

**Synonyms:** MI0021881; MIMAT0025096; MIR6353 mmu-mir-6353; mmu-miR-6353

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)