

Product datasheet for **MK300831**

MIR6396 Mouse qPCR Template Standard (MI0021931)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: GGGGCCTTAGCTCCTAGTCC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025147

Synonyms: MI0021931; MIMAT0025147; MIR6396 mmu-mir-6396; mmu-miR-6396

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)