

## Product datasheet for **MK300872**

### **MIR669a-8 Mouse qPCR Template Standard (MI0014066)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Forward Sequence:</b>	AGTTGTGTGTCATGTTTCATG
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0003477
<b>Synonyms:</b>	MI0014066; MIMAT0003477; MIR669a-8 mmu-mir-669a-8; mmu-miR-669a-5p
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
<b>Quality Control:</b>	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)