

Product datasheet for **MK300915**

MIR6931 Mouse qPCR Template Standard (MI0022778)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGGGTGGGGAGTGGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027762
Synonyms:	MI0022778; MIMAT0027762; MIR6931 mmu-mir-6931; mmu-miR-6931-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)