

## Product datasheet for **MK300977**

### **MIR6992 Mouse qPCR Template Standard (MI0022840)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** TGCCTGTGATGGTTTGGCTG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0027886

**Synonyms:** MI0022840; MIMAT0027886; MIR6992 mmu-mir-6992; mmu-miR-6992-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)