

Product datasheet for **MK301016**

MIR7031 Mouse qPCR Template Standard (MI0022880)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CCTGAGAGGCCTGAAGGGTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0027966

Synonyms: MI0022880; MIMAT0027966; MIR7031 mmu-mir-7031; mmu-miR-7031-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)