

## Product datasheet for **MK301035**

### **MIR7049 Mouse qPCR Template Standard (MI0022898)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** TCTGAAGGCAGTCAGGTAGGC

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0028002

**Synonyms:** MI0022898; MIMAT0028002; MIR7049 mmu-mir-7049; mmu-miR-7049-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)