

Product datasheet for **MK301081**

MIR7094-2 Mouse qPCR Template Standard (MI0022945)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCTGGAGAGGATACAGGGTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0028096
Synonyms:	MI0022945; MIMAT0028096; MIR7094-2 mmu-mir-7094-2; mmu-miR-7094b-2-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)