

## Product datasheet for **MK301082**

### MIR7115 Mouse qPCR Template Standard (MI0022966)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** TCGGGGGTTGTGGTGCCGA

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0028127

**Synonyms:** MI0022966; MIMAT0028127; MIR7115 mmu-mir-7115; mmu-miR-7115-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)