

## Product datasheet for **MK301088**

### **MIR7211 Mouse qPCR Template Standard (MI0023706)**

#### **Product data:**

|                   |                                                                                                                                                                                                       |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Template Standards                                                                                                                                                                               |
| Forward Sequence: | TCTTCCCTCTGCCACTCC                                                                                                                                                                                    |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                   |
| ACCN:             | MIMAT0028390                                                                                                                                                                                          |
| Synonyms:         | MI0023706; MIMAT0028390; MIR7211 mmu-mir-7211; mmu-miR-7211-5p                                                                                                                                        |
| Components:       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control:  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| Storage:          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |



[View online »](#)