

## Product datasheet for **MK301095**

### MIR7218 Mouse qPCR Template Standard (MI0023713)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGCAGGGTTTAGTGTAGAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0028404
Synonyms:	MI0023713; MIMAT0028404; MIR7218 mmu-mir-7218; mmu-miR-7218-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)