

Product datasheet for **MK301100**

MIR7223 Mouse qPCR Template Standard (MI0023718)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGCTTTGCTGAGGTGTCG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0028414

Synonyms: MI0023718; MIMAT0028414; MIR7223 mmu-mir-7223; mmu-miR-7223-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)