

Product datasheet for **MK301105**

MIR7228 Mouse qPCR Template Standard (MI0023723)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGCGACCTGAACAGATG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0028424
Synonyms:	MI0023723; MIMAT0028424; MIR7228 mmu-mir-7228; mmu-miR-7228-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)