

Product datasheet for **MK301121**

MIR7578 Mouse qPCR Template Standard (MI0024714)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CATGGCTCTGTCTTCTGCCT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0029578

Synonyms: MI0024714; MIMAT0029578; MIR7578 mmu-mir-7578; mmu-miR-7578

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)