

Product datasheet for **MK301128**

MIR7652 Mouse qPCR Template Standard (MI0024992)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TAAGGGCACAGGGATTTTCAG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0029810

Synonyms: MI0024992; MIMAT0029810; MIR7652 mmu-mir-7652; mmu-miR-7652-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)