

## Product datasheet for **MK301140**

### **MIR7664 Mouse qPCR Template Standard (MI0025004)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** CAGTTTACCTGTCTATGAG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0029834

**Synonyms:** MI0025004; MIMAT0029834; MIR7664 mmu-mir-7664; mmu-miR-7664-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)