

Product datasheet for **MK301153**

MIR7676-1 Mouse qPCR Template Standard (MI0025017)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGGGCAGTATGATGGCCTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0029860
Synonyms:	MI0025017; MIMAT0029860; MIR7676-1 mmu-mir-7676-1; mmu-miR-7676-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)