

## Product datasheet for **MK301189**

### **MIR8110 Mouse qPCR Template Standard (MI0026042)**

#### **Product data:**

|                   |                                                                                                                                                                                                       |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Template Standards                                                                                                                                                                               |
| Forward Sequence: | AAGCGTGGATTGGGGGGG                                                                                                                                                                                    |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                   |
| ACCN:             | MIMAT0031416                                                                                                                                                                                          |
| Synonyms:         | MI0026042; MIMAT0031416; MIR8110 mmu-mir-8110; mmu-miR-8110                                                                                                                                           |
| Components:       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control:  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| Storage:          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |



[View online »](#)