


Product datasheet for **MO18004-100**

CD34 Mouse Monoclonal Antibody [Clone ID: ICO115]

Product data:

Product Type:	Primary Antibodies
Clone Name:	ICO115
Applications:	IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	CD34 molecule
Database Link:	Entrez Gene 947 Human
Synonyms:	Hematopoietic progenitor cell marker


Product images:



CD34 staining of paraffin-embedded human capillary hemangioma.



[View online »](#)



CD34 staining of paraffin-embedded human capillary hemangioma.