

Product datasheet for **MO20000-100**

Synaptophysin (SYP) Mouse Monoclonal Antibody [Clone ID: 27G12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	27G12
Applications:	IHC
Reactivity:	Human, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Supernatant
Gene Name:	Homo sapiens synaptophysin (SYP)
Database Link:	Entrez Gene 24804 Rat Entrez Gene 6855 Human
Synonyms:	Neuroendocrine Marker
Protein Families:	Secreted Protein, Transmembrane

Product images:



[View online »](#)

