


Product datasheet for **MO20002-100**

PGP9.5 (UCHL1) Mouse Monoclonal Antibody [Clone ID: 10A1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	10A1
Applications:	IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Formulation:	State: Supernatant
Gene Name:	ubiquitin C-terminal hydrolase L1
Database Link:	Entrez Gene 7345 Human
Synonyms:	UCH-L1, PGP 9.5, Ubiquitin thioesterase L1, Neuron cytoplasmic protein 9.5

Product images:



[View online »](#)

